

**FREE** **DOM** BathTubs

FreedomBathtubs.com

ACCESSIBILITY PROFESSIONALS **AP**

**AQUARITE**

**Installation Manual**


**APWHW30**




# Table of Contents

Diagrams	.....	3
Notice to Installers	.....	4
Tools Needed for Installation	.....	4
Tub Installation & Assembly	.....	4
Drain Installation & Assembly	.....	6


Left Hand Tub




Right Hand Tub



Left Hand Drain Tub



Right Hand Drain Tub



- Height to step in: 9.5"
- Inside height of tub: 18.5"
- Overflow depth: 16" to drain
- Tub height might vary +-1" due to adjustable leveling feet.

## **NOTICE TO INSTALLERS:**

- ✓ Tub needs to be leveled in its final resting position, mark the placement, and then level. All feet must be touching the floor.
- ✓ It is advised to have tub onsite before preparing drain location.
- ✓ Make sure the floor can withstand the weight of the product.
- ✓ This is a free standing tub. You must account for the length of the skirt if tiling in the tub.


*Manufacturer reserves the right to change specs or features at anytime. Please check to confirm details. 1-866-783-2661*

## **TOOLS NEEDED FOR INSTALLATION:**




## **INSTALLATION**

1. Walls and floor must be clear of any debris. Note: Tub comes with 90 degree drain; make sure your drain location will work with the drain.




**\*Please allow room for the skirt after installation.**


2. Remove the skirt off the tub for easy access.




- 3 Place the tub in its final resting place; adjust all the legs to be touching the floor. Level the tub in place.



2. With the level, check the outside of the tub to be leveled. Adjust the feet for leveling.



3. Fix the bathtub; tighten the nuts on the leveling feet.



4. Water test:  
After installation, the tub has to be checked for leaks. Close and lock the glass door. Fill the tub to the overflow control.

**NOTE:** Any gaps between the skirt and the tub can be caulked with latex caulk. The caulk can be removed and replaced when the tub needs service.

# DRAIN ASSEMBLY INFORMATION

**Note:** One 8in section of Schedule 40, 1 ½in drain pipe for item #8 of water bath waste system will be needed in addition to the drain kit included with the tub.

**Step 1:** Cut the 1 1/2" Schedule 40 pipes (Item 8 and Item 10) to the desired lengths to fit the tub.

**Step 2:** Using the proper solvent cement, assemble the bath drain pipes and fittings (Item 7 thru Item 11).

**Step 3:** Place foam gasket (Item 12) onto overflow elbow (Item 11).

**Step 4:** Insert overflow elbow and foam gasket assembly into the tub overflow hole from the back of the tub.

**Step 5:** Screw on retainer nut (Item 13) from inside of tub. hand-tighten only.

**Step 6:** Secure drain elbow (Item 7) and drain flange gasket (Item 6) under the tub by screwing in the strainer body (Item 5) through tub drain hole from inside the tub. If not testing the bath waste, skip to Step 8.


**Step 7:** Screw the strainer test plug (Item 15) into the strainer body (Item 5) and test the system. Once testing is complete, remove the strainer test plug (Item 15).

**Step 8:** Cut away the orange test membrane from the overflow elbow (Item 11).

**Step 9:** Snap on the overflow plate (Item 14) with the slot facing downward.

**Step 10:** Install the stopper assembly (Item 1 thru Item 4) into the strainer body (Item 5).

**Note:** Items 1 thru 4 may vary depending on the type of stopper assembly being used.



## Accessibility Professionals Inc.

4 Alliance Blvd. Unit 123

Barrie, ON L4M 7G3

Toll Free: [877-947-7769](tel:877-947-7769)

Toll Free Fax: [877-794-7850](tel:877-794-7850)

[www.FreedomBathtubs.com](http://www.FreedomBathtubs.com)